



Cost-effective,  
 needleless system  
 with sterile  
 connection and  
 single donor exposure

## Neonatal/Pediatric Blood Administration Set

The **Marian Medical Blood Administration Set** is a syringe set designed for the preparation and delivery of small aliquots of whole blood, red cells and platelets with a 150-micron filter. It is a cost effective, needleless system that maximizes the usage of adult size units of blood with sterile connection and single donor exposure for the neonatal/pediatric patient.

- + *Non-DEHP, Latex-free and BPA-free blood delivery system*
- + *Sterile, pyrogen-free blood path*
- + *Compatible with sterile connecting technology*
- + *Provides a closed blood delivery system*
- + *Limits the multiple donor exposure to the neonate/pediatric patient*
- + *150 micron filter in accordance with AABB blood administration requirements*
- + *Needleless system with protected spike for safe handling*



### Neonatal/Pediatric Syringe Set

The system includes a BDTM syringe that is accurately identified and compatible with all IV pumps..

	Syringe	Micron Filter	Qty/Box
<b>BASN-30/150</b>	30 cc	150	25
<b>BASN-60/150</b>	60 cc	150	25

\*BD is a registered trademark of Becton Dickinson and Co.



**VITAL CARE PRODUCTS, INC.**

3046 Brecksville Road + Suite 2 + Richfield, Ohio 44286  
 1 800 837 5511 + [www.vitalcareproducts.com](http://www.vitalcareproducts.com)

*critical care solutions delivered full circle*

