



Cost-effective,
 needleless system
 with sterile
 connection and
 single donor exposure

Neonatal Exchange Transfusion Tray

The **Marian Medical Exchange Transfusion Tray** is a valuable tool in standardizing protocols for quality patient care. It is a cost effective, needleless system that maximizes the use of sterile connection and single donor exposure for the neonatal patient.

Neonatal Exchange Transfusion Tray includes:

- 1 CSR Wrap
- 4 Documentation Sheets
- 1 Blood Tubing 67" w/bag spike and 150 micron filter
- 1 Blood Waste Bag
- 1 Hanging Hook for blood waste bag
- 1 Extension Set 33"L with MLL and FLL on opposite end
- 1 Shuttle Set
- 1 Syringe 10 ml
- 1 Syringe 30 ml



MR-ETT-NW



Neonatal Exchange Transfusion Tray

MR-ETT-NW Single use only Latex & DEHP Free 1 box

** Directions are supplied to help explain the use of our components.
 Procedural details may vary according to clinical judgment as well as hospital.*




VITAL CARE PRODUCTS, INC.
 3046 Brecksville Road + Suite 2 + Richfield, Ohio 44286
 1 800 837 5511 + www.vitalcareproducts.com

critical care solutions delivered full circle




Neonatal Exchange Transfusion Tray

Single Catheter Setup with Shuttle Set



* MLL port and FLL Connection can be disconnected for use with **dual** catheter setup.

Single Catheter Connection Completed



Warning: Umbilical catheters should NOT be left in place between procedures.

Caution: Blood administration set should be used if protectors are not in place. Do not allow air to be trapped in the line/set. Do not puncture blood chamber or tubing because that could cause an air embolism.