

**The Leader In  
PFT Filters!**


# PulmoGUARD™

*Fits most PFT devices without adapters*

During a Pulmonary Function Test, bacteria and viruses can be expelled, posing serious risk to patients and health care workers alike. With a cross-contamination efficiency of 99.999%, PulmoGuard is safe, highly efficient...and very affordable.

## PulmoGUARD™

INFECTION PROTECTION



- Meets ATS air resistance standards
- Dead space of 60cc
- Pleated media for better functionality
- Compact, non-intimidating design
- Transparent for easy, 100% visual inspection
- Fits most spirometers with no adapter

*Please see back for Technical Specifications and uses.*


### PulmoGuard Part Numbers & Packaging

Part No.	Description	Qty	Packaging
29-7950-050	PulmoGuard Filter	50	Box
29-7950-100	PulmoGuard Filter	100	Box
29-7950-200	PulmoGuard Filter	200	4 Boxes/50
29-7950-300	PulmoGuard Filter	300	Box

*Please see back for Kit options*

## PulmoGUARD II™

INFECTION PROTECTION



- All the benefits of the original PulmoGuard  
**PLUS –**
- Fits many PFT instruments with no adapter
  - Dead space of 64cc

*Please see back for Technical Specifications and uses.*

### PulmoGuard II Part Numbers & Packaging

Part No.	Description	Qty	Packaging
29-7920-050	PulmoGuard II Filter	50	Box
29-7920-100	PulmoGuard II Filter	100	Box
29-7920-200	PulmoGuard II Filter	200	4 Boxes/50
29-7920-300	PulmoGuard II Filter	300	Box

*Please see back for Kit options*

Your Dealer is:

**VITAL CARE PRODUCTS, INC.**  
 sales@vitalcareproducts.com  
 800 837 5511

**sdi**

DIAGNOSTICS

*The Spirometry Source*

## PULMOGUARD TECHNICAL SPECIFICATIONS

PulmoGuard fits the following instruments without an adapter:  
Adapters available for instruments not listed.

**Assess** - Peak Flow Meter  
**Burdick** - Spirosense  
**Circadian** - Models (if using black Tamarac pneumotach, need 7975)  
**Collins** - Survey, Survey II, Survey III, Survey PLUS, Eagle, Eagle II  
**Cybermedics** - Spinnaker TL, SL, CM3, CM310 (if white or black cone use 7980)  
**Cybermedics** - Moose  
**Godart** - All Models  
**Gould** - 5004  
**Healthscan** - Assess Peak Flow Meter  
**Jaeger** - All Models except Spiro Pro  
**Jones** - All Models except Satellite (small end of filter)  
**Koko** - All Models  
**Med Graphics** - 1070, 1085 (w/Prevent, use PG II; if no Prevent, 7970-spirometry, 7961-vol/diff)  
**Med Science** - 570, 6257  
**Med Science** - 3000, 4000 (spirometry-no adapt./vol-7961/diff-7960B)  
**Med Systems** - Spirolite 101 & 201 (small end of filter), 343  
**New Med** - 7000 I (Spirobyte)  
**PK Morgan** - Benchmark (spirometry-7978, vol&diff-7960B), Flexi Flow (small end of filter), Trans Flow  
**Pneumedics** - All Models (small end of filter)  
**Riko** - AS300  
**S&M** - CF-150B, GS520, Keystone III, Portascreen, Pulmoscreen  
**ViaSys/Sensormedics** - Horizon, Ohio, 922, 2130, 2200, 2400, 2450, 2800, 6200, 2100 (spirometry-7961, vol-7960B, diff-7960B), Vmax (can use PGII w/7976 or 7977)  
**Spirolink** - All Models  
**Spirometrics** - LTE, 2500 LTE  
**Spiroscreen** - All Models (small end of filter)  
**Spirotech** - All Models  
**Tamarac** - All Models except Presto Flash

## PULMOGUARD II TECHNICAL SPECIFICATIONS

PulmoGuard II fits the following instruments without an adapter:  
Adapters available for instruments not listed.

**Breon** - All Models  
**CDX** - Spiro 110 (remove black fitting), Spiro 850  
**Cavitron/Donti** - All Models  
**Clement-Clarke** - Mini-Wright Peak Flow Meter  
**Collins** - All DS series pulmonary labs with UBV (use 7982 adapter for hose), CPL, GS, SurveyTach, Eaglet, Hawk  
**Cosmed** - All Models  
**Creative Biomedics** - All Models  
**Futuremed** - Discovery, Spirobank, Spirovision III  
**Medical Equip Designs** - All Models (remove black fitting)  
**Med Graphics** - preVent, PF/DX (needs to be reprogrammed if no filter before)  
**Micro Medical/Direct** - All Models  
**MIR** - All Models  
**Multispiro** - All Models except PC (remove black fitting)  
**nSpiro** - Wright Peak Flow Meter  
**Ohio** - All Models  
**PK Morgan** - 232  
**Puritan-Bennett** - PS600, 900, 900L, VS400  
**Riko** - All Models except AS300  
**S&M** - Keystone  
**SDI** - SBOffice, SBG, spirotel, Spirolab, Spirolab II, Astra series with 7983 adapter  
**Smiths Medical** - All models  
**Spirometrics** - All Models except LTE, 2500 LTE  
**VacuMed** - Discovery, Micro, MicroPlus, MicroGP/DL, MicroLoop, Pony, Spiroscreen, Wizard  
**ViaSys** - MicroLoop, MicroLab  
**Vitalograph** - All Models except Presto Flash  
**Wright Professional Peak Flow Meter**

## PulmoGuard Kits

Part No.	Description	Qty	Packaging
29-7951-050	PulmoGuard Filter + Cardboard Mouthpiece	50	Box
29-7951-100	PulmoGuard Filter + Cardboard Mouthpiece	100	Box
29-7953-040	PulmoGuard Filter + Cardboard Mouthpiece + The Klip	40	Box
29-7953-080	PulmoGuard Filter + Cardboard Mouthpiece + The Klip	80	Box
29-7954-050	PulmoGuard Filter + The Klip	50	Box
29-7954-100	PulmoGuard Filter + The Klip	100	Box
29-7957-040	PulmoGuard Filter + Comfit-D Disposable Mouthpiece	40	Box
29-7957-080	PulmoGuard Filter + Comfit-D Disposable Mouthpiece	80	Box
29-7958-040	PulmoGuard Filter + Comfit-D Disposable Mouthpiece + The Klip	40	Box
29-7958-080	PulmoGuard Filter + Comfit-D Disposable Mouthpiece + The Klip	80	Box

## PulmoGuard II Kits

Part No.	Description	Qty	Packaging
29-7921-050	PulmoGuard II Filter + 1" O.D. Mouthpiece	50	Box
29-7921-100	PulmoGuard II Filter + 1" O.D. Mouthpiece	100	Box
29-7922-050	PulmoGuard II Filter + 1-1/16" O.D. Mouthpiece	50	Box
29-7922-100	PulmoGuard II Filter + 1-1/16" O.D. Mouthpiece	100	Box
29-7923-040	PulmoGuard II Filter + 1" O.D. Mouthpiece + The Klip	40	Box
29-7923-080	PulmoGuard II Filter + 1" O.D. Mouthpiece + The Klip	80	Box
29-7924-050	PulmoGuard II Filter + The Klip	50	Box
29-7924-100	PulmoGuard II Filter + The Klip	100	Box
29-7925-040	PulmoGuard II Filter + 1-1/16" O.D. Cardboard Mouthpiece + The Klip	40	Box
29-7925-080	PulmoGuard II Filter + 1-1/16" O.D. Cardboard Mouthpiece + The Klip	80	Box
29-7927-040	PulmoGuard II Filter + Comfit-D Disposable Mouthpiece	40	Box
29-7927-080	PulmoGuard II Filter + Comfit-D Disposable Mouthpiece	80	Box
29-7928-040	PulmoGuard Filter + Comfit-D Disposable Mouthpiece + The Klip	40	Box
29-7928-080	PulmoGuard Filter + Comfit-D Disposable Mouthpiece + The Klip	80	Box